

	Product #	Description	Storage Temperature	Fill Size	Price
<b>Animal Serum</b>	FB-01	Fetal Bovine Serum	-20°C	500mL	Prices available upon request
	FB-02	Fetal Bovine Serum, Heat Inactivated	-20°C	500mL 100mL	Prices available upon request
	FB-03	Fetal Bovine Serum, Dialyzed	-20°C	500mL 100mL	Prices available upon request
	FB-04	Fetal Bovine Serum, Charcoal / Dextran Treated	-20°C	500mL 100mL	Prices available upon request
	FB-05	Fetal Bovine Serum, ES Cell Tested	-20°C	500mL 100mL	Prices available upon request
	FB-06	Fetal Bovine Serum, Low IgG	-20°C	500mL	Prices available upon request
	NC-03	Newborn Calf Serum	-20°C	500mL 100mL	\$45.00 \$15.00
	BC-04	Bovine Calf Serum	-20°C	500mL 100mL	\$40.00 \$10.00
	IS-06	Bovine Calf Serum, Iron Supplemented	-20°C	500mL 100mL	\$50.00 \$15.00
	DH-05	Donor Horse Serum	-20°C	500mL 100mL	\$42.00 \$15.00
	DH-05A	Agamma Horse Serum	-20°C	500mL	\$140.00
	DH-07	Donor Horse Serum, Heat Inactivated	-20°C	500mL 100mL	\$42.00 \$15.00
	MS-07	Mouse Serum	-20°C	100mL	\$250.00
	NG-11	Goat Serum	-20°C	500mL 100mL	\$100.00 \$25.00

Please call for quotes on 12 bottles or more.

<b>Serum Treatments</b>	Omega Scientific, Inc. offers a range of serum treatments. For most, these are available as stock items for many types of animal sera. Others may be ordered on request.				
	<b>Heat Inactivation</b>	Serum is heated to and maintained at a control temperature of 56 degrees Celsius for 30 minutes.			
	<b>Dialyzed</b>	Serum is Dialyzed against physiological saline with a 10,000 molecular weight cutoff. The glucose level is monitored until it is less than 5.0 mg/dL. Exhaustive dialysis is not performed in order to prevent precipitation and inactivation of peptides.			
	<b>Gamma Irradiated</b>	Serum is irradiated with a minimum dose of 2.5 megarads.			
	<b>Reduced IgG</b>	Fetal Bovine Serum is processed by a proprietary chromatography method until the IgG level is less than 5µg/mL. The method does not remove any of the other serum proteins or growth factors.			

	Product #	Description	Storage Temperature	Fill Size	Price
<b>Human Serum</b>	HS-20	Human AB Serum, Male Donors	-20°C	100mL	\$95.00
	HS-30	Human Serum, Male Donors	-20°C	100mL	\$85.00
	<b>Human Serum Testing</b>	Human Serum is collected from companies whose screening practices are believed to conform to all applicable FDA regulations. Specifically we are informed that each individual donor unit before pooling is tested and found negative for Hepatitis B (HBsAG), Hepatitis C (Anti-HCV) and Human Immunodeficiency Virus (HIV) antibodies. Final product is tested and found negative for (HBsAG), Hepatitis A and (HIV). Nevertheless, products of human origin should be considered potentially infectious and handled accordingly.			

	Product #	Description	Storage Temperature	Fill Size	Price
<b>Human Albumin</b>	AB-05	Human Albumin, 5% Solution, Therapeutic Grade	Room Temperature	50mL	Prices available upon request
				100mL	
				250mL	
				500mL	
AB-20	Human Albumin, 20% Solution, Therapeutic Grade	Room Temperature	50mL	Prices available upon request	
			100mL		
AB-25	Human Albumin, 25% Solution, Therapeutic Grade	Room Temperature	50mL	Prices available upon request	
			100mL		
	<b>Human Albumin Testing</b>	Human Albumin has been manufactured starting from plasma units coming from the United States. These units have been analyzed and found non-reactive through the detection tests for Hepatitis B surface antigen (HBsAg), antibodies against human immunodeficiency viruses 1 and 2, HIV 1-2, and antibodies against Hepatitis C virus (HCV), using FDA licensed tests. Moreover, the pool of source plasma of fractionation from which the products come from has been analyzed and found non-reactive through the same techniques mentioned above.			

	Product #	Description	Storage Temperature	Fill Size	Price
<b>Bovine Albumin</b>	BA-90	Bovine Albumin	Room Temperature	100 grams 500 grams 1 kg	Prices available upon request

	Product #	Description	Storage Temperature	Fill Size	Price
<b>Tissue Culture Media</b>	RP-18	RPMI 1640 Medium (1X) w/o L-Glutamine with 25mM HEPES	2-8°C	500mL	\$17.50
				1000mL	\$35.00
	RP-19	RPMI 1640 Medium (1X) with L-Glutamine	2-8°C	500mL	\$8.95
				1000mL	\$15.95
	RP-20	RPMI 1640 Medium (1X) w/o L-Glutamine	2-8°C	500mL	\$8.95
				1000mL	\$15.95
	DM-21	DMEM Low Glucose (1000mg/L) w/o L-Glutamine	2-8°C	500mL	\$8.95
				1000mL	\$15.95
	DM-22	DMEM High Glucose (4500mg/L) w/o L-Glutamine	2-8°C	500mL	\$8.95
				1000mL	\$15.95
DM-25	DMEM Ham's F-12 (1:1 ratio) w/o L-Glutamine	2-8°C	500mL	\$19.95	
ME-23	MEM Earle's Salts (1X) w/o L-Glutamine	2-8°C	500mL	\$8.95	
	PR-10	PRBS (10X) w/o Calcium & Magnesium	2-8°C	500ml	\$12.00

	PB-10	PBS (10X) w/o Calcium & Magnesium	2-8°C	500mL	\$12.00
	PB-24	PBS (1X) w/o Calcium & Magnesium	2-8°C	500mL	\$10.00
				1000mL	\$20.00
	IM-25	IMDM (1X) 25mM HEPES w/o L-Glutamine	2-8°C	500mL	\$19.00
				1000mL	\$38.00
	HB-10	HBSS w/o Calcium & Magnesium	2-8°C	1000mL	\$15.00
	<b>Product #</b>	<b>Description</b>	<b>Storage Temperature</b>	<b>Fill Size</b>	<b>Price</b>
<b>Liquid Antibiotics &amp; Supplements</b>	PS-20	Penicillin Streptomycin Solution 100X 10,000 units Penicillin-G/mL 10,000 mcg Streptomycin/mL in normal saline	-20°C	100mL	\$9.95
	PG-30	L-Glutamine Pen-Strep Solution 100X L-Glutamine 200mM 10,000 units Penicillin-G/mL 10,000 mcg Streptomycin/mL in normal saline	-20°C	100mL	10.95
	AA-40	Antibiotic Antimycotic Solution 100X 10,000 units Penicillin-G/mL 10,000 mcg Streptomycin/mL 25mcg/mL Fungizone™ in normal saline	-20°C	100mL	\$13.95
	GS-60	L-Glutamine Solution, 200mM	-20°C	100mL	\$9.95
	GT-10	Gentamicin Sulfate Solution 10 mg/mL in water	Room Temperature	6x25mL	\$55.00 (per box)
	GT-50	Gentamicin Sulfate Solution 50 mg/mL in water	Room Temperature	6x25mL	\$175.00 (per box)
	FG-70	Fungizone™	-20°C	100mL	\$14.95
	HG-80	Hygromycin B	Room Temperature	20mL	\$95.00
	NE-01	Nonessential Amino Acids 10mM NEAA (100X)	Room Temperature	100mL	\$9.95
	SP-90	Sodium Pyruvate 100mM	2-8°C	100mL	\$8.95
	HB-20	Hepes Buffer 1M	2-8°C	100mL	\$15.50
	TE-91	Trypsin EDTA 1X 0.5 g/L Trypsin (1:250) 0.2 g/L EDTA in HBSS	-20°C	100mL	\$6.95
	™Trademark of E.R. Squibb				
	<b>Product #</b>	<b>Description</b>	<b>Storage Temperature</b>	<b>Fill Size</b>	<b>Price</b>
<b>Powdered Antibiotics</b>	CB-01	Carbenicillin	2-8°C	5 grams	\$95.00
	CF-02	Cefotaxime	2-8°C	5 grams	\$195.00
	GT-03	Gentamicin	Room Temperature	5 grams	\$75.00
	GN-04	G-418 Sulfate	Room Temperature	5 grams	\$225.00
1 gram				\$55.00	
	<b>Product #</b>	<b>Description</b>			<b>Price</b>
<b>Aerosol Barrier Tips</b>	1502	Aerosol Tips 0.5-10µL working volume (For Eppendorf Pipettors) Sterile, Racked 960/pack			\$65.00/Pack
	1508	Aerosol Tips 0.1-10µL working volume Sterile, Racked, Universal Fit 960/pack			\$65.00/Pack
	1421	Aerosol Tips 1-30µL working volume Sterile, Racked, Universal Fit 960/pack			\$65.00/Pack
	1422	Aerosol Tips 5-200µL working volume Sterile, Racked, Universal Fit 960/pack			\$65.00/Pack
	1491	Aerosol Tips 5-250µL working volume Sterile, Racked, Universal Fit			\$65.00/Pack

		960/pack	
	1420	Aerosol Tips 200-1000 $\mu$ L working volume Sterile, Racked, Universal Fit 1000/pack	\$65.00/Pack
	Product #	Description	Price
Gel Loading Tips	1376	Flat Gel Tip 0.37mm OD Natural, Flat End Tips, Racked 200/rack, 4 racks per case	\$115.00/Case
	1735	Flat Gel Tip 0.1-10 $\mu$ L 0.17mm OD Natural, Flat End Tips, Racked 200/rack, 4 racks per case	\$159.00/Case
	1379	Round Gel Tip 0.57mm OD Natural, Racked 200/rack, 4 racks per case	\$69.00/Case
	Product #	Description	Price
Liquid Handling Plastics	1505	Natural Tips 0.1-10 $\mu$ L, Bulk Universal Fit, Beveled 1000/pack, 10 packs per case	\$129.00/Case
	1503	Natural Tips 0.1-10 $\mu$ L, Racked Universal Fit, Beveled 960/pack, 5 packs per case	\$139.00/Case
	1504	Natural Tips 0.1-10 $\mu$ L, Racked Sterile, Universal Fit, Beveled 960/pack, 5 packs per case	\$149.00/Case
	2221	Natural Tips 0.5-10 $\mu$ L, Racked Sterile, For Eppendorf, Beveled 1000/pack, 5 packs per case	\$155.00/Case
	2222	Natural Tips 0.5-10 $\mu$ L, Racked For Eppendorf, Beveled 1000/pack, 5 packs per case	\$145.00/Case
	1531	Yellow Tips 1-200 $\mu$ L, Bulk Universal Fit, Beveled 1000/pack, 10 packs per case	\$105.00/Case
	1532	Yellow Tips 1-200 $\mu$ L, Racked Universal Fit, Beveled 960/pack, 5 packs per case	\$115.00/Case
	1533	Yellow Tips 1-200 $\mu$ L, Racked Sterile, Universal Fit, Beveled 960/pack, 5 packs per case	\$125.00/Case
	UY0200-MS	Yellow Tips 1-200 $\mu$ L, Beveled Universal Fit, Stacked Racked (5x960)	\$89.00/Case
	UY0200-MSS	Yellow Tips 1-200 $\mu$ L, Beveled Sterile, Universal Fit, Stacked Racked (5x960)	\$105.00/Case
	3003	Yellow Tips 1-200 $\mu$ L, Starter Kit Quickshot Reloading System Universal Fit, Beveled 1 magazine & 10 empty racks per kit	\$49.00/Kit
	3025	Yellow Tips 1-200 $\mu$ L, Magazine Quickshot Reloading System Universal Fit, Beveled 960 tips/magazine, 10 magazines per case	\$210.00/Case

900 tips/magazine, 10 magazines per case

	3046	Yellow Tips 1-200 $\mu$ L, Magazine Quickrack Reloading System Universal Fit, Beveled 960 tips/magazine, 10 magazines per case	\$90.00/Case
	1004	Blue Tips 200-1000 $\mu$ L, Bulk Universal Fit 1000/pack, 10 packs per case	\$120.00/Case
	1008	Blue Tips 200-1000 $\mu$ L, Racked Universal Fit 1000/pack, 5 packs per case	\$129.00/Case
	1013	Blue Tips 200-1000 $\mu$ L, Racked Sterile, Universal Fit 1000/pack, 5 packs per case	\$145.00/Case
	<b>Product #</b>	<b>Description</b>	<b>Price</b>
CELLine Cell Culture Flasks	90005	CELLine CL 350 Flasks Disposable Two-Compartment Bioreactor for Suspension Cells 350 mL media volume, 5 mL culture volume 3 flasks per case	\$535.00/Case
	90010	CELLine CL 1000 Flasks Disposable Two-Compartment Bioreactor for Suspension Cells, 1000 mL media volume, 15 mL culture volume 5 flasks per case	\$739.00/Case
	90025	CELLine AD 1000 Flasks Disposable Two-Compartment Bioreactor with Matrix Inlay for Anchorage-Dependent Cells 1000 mL media volume, 15 mL culture volume 3 flasks per case	\$536.00/Case
	<b>Product #</b>	<b>Description</b>	<b>Price</b>
Liquid Handling / Manual Pipetting	155 000	Pipetboy, Classic	\$325.00
	155 015	Pipetboy, Green	\$325.00
	155 016	Pipetboy, Red	\$325.00
	155 017	Pipetboy, Blue	\$325.00
	155 018	Pipetboy, Transparent	\$325.00
	155 019	Pipetboy, Purple	\$325.00
	153 015	Filter Sterile, 0.45 $\mu$ m 10 per pack	\$47.00
	153 005	9V NiMH Battery	\$30.00
	<b>Product #</b>	<b>Description</b>	<b>Price</b>
Microcentrifuge Tubes	1151	Microcentrifuge Tube, 1.7mL Graduated, Flat Top Natural, Boilproof 500/pack, 10 packs per case	\$109.00/Case
	1115	Microcentrifuge Tube, 0.65mL Graduated, Flat Top Natural, Boilproof 1000/pack, 10 packs per case	\$215.00/Case
	1121	Microcentrifuge Tube, 0.65mL Graduated, Flat Top Rainbow, Boilproof 1000/pack, 10 packs per case	\$220.00/Case
	1159	Microcentrifuge Tube, 1.7mL Graduated, Flat Top Rainbow, Boilproof 500/pack, 10 packs per case	\$119.00/Case
	<b>Product #</b>	<b>Description</b>	<b>Price</b>
	1516	Ultra PCR Tubes, 0.65mL Thin-Wall, Natural RNase/DNase Free 1000/pack, 10 packs per case	\$299.00/Case

<b>Ultra PCR Products</b>	1695	Ultra PCR Tubes, 0.2mL Thin-Wall, Natural RNase/DNase Free 1000/pack, 10 packs per case	\$299.00/Case
	1633	Ultra Strip PCR Tubes, 0.2mL Thin-Wall, Natural, with caps RNase/DNase Free 8 tubes/strip, 5 bags of 24 strips per case	\$89.00/Case
	<b>Product #</b>	<b>Description</b>	<b>Price</b>
<b>Microflex Gloves</b>	EV2050	Evolution One Latex, Powder Free 100/box, 10 boxes per case Available in XS, S, M, L, XL	\$89.00/Case
	MF300	Diamond Grip Latex, Powder Free 100/box, 10 boxes per case Available in XS, S, M, L, XL	\$89.00/Case
	DGP350	Diamond Grip Plus Latex, Powder Free 100/box, 10 boxes per case Available in XS, S, M, L, XL	\$80.00/Case
	PG199	Power Grip Latex, Low Powder 100/box, 10 boxes per case Available in XS, S, M, L, XL	\$65.00/Case
	CT144	Color Touch Latex, Low Powder 100/box, 10 boxes per case Available in S, M, L (coded by color)	\$69.00/Case
	UL315	Ultra One Latex, Powder Free 50/box, 10 boxes per case Available in S, M, L, XL	\$105.00/Case
	SY911	Synetron Latex, Powder Free 50/box, 10 boxes per case Available in S, M, L, XL	\$105.00/Case
	NP144	Neopro Chloroprene, Powder Free 100/box, 10 boxes per case Available in XS, S, M, L, XL	\$100.00/Case
	NO123	Nitron One Nitrile, Low Powder 100/box, 10 boxes per case Available in S, M, L, XL	\$80.00/Case
	SP690	Supreno Nitrile, Powder Free 100/box, 10 boxes per case Available in XS, S, M, L, XL	\$100.00/Case
	Please designate size with catalog number when ordering.		
	<b>Product #</b>	<b>Description</b>	<b>Price</b>
	MDS195013	Aloetouch Latex, Powder Free 100/box, 10 boxes per case Extra-Small	\$68.00/Case
	MDS195014	Aloetouch Latex, Powder Free 100/box, 10 boxes per case Small	\$68.00/Case
	MDS195015	Aloetouch Latex, Powder Free 100/box, 10 boxes per case Medium	\$68.00/Case
	MDS195016	Aloetouch Latex, Powder Free 100/box, 10 boxes per case Large	\$68.00/Case

		Large	
<b>Aloetouch Gloves</b>	MDS195017	Aloetouch Latex, Powder Free 100/box, 10 boxes per case Extra-Large	\$68.00/Case
	MDS195083	Aloetouch Nitrile, Powder Free 100/box, 10 boxes per case Extra-Small	\$89.50/Case
	MDS195084	Aloetouch Nitrile, Powder Free 100/box, 10 boxes per case Small	\$89.50/Case
	MDS195085	Aloetouch Nitrile, Powder Free 100/box, 10 boxes per case Medium	\$89.50/Case
	MDS195086	Aloetouch Nitrile, Powder Free 100/box, 10 boxes per case Large	\$89.50/Case
	MDS195087	Aloetouch Nitrile, Powder Free 100/box, 10 boxes per case Extra-Large	\$89.50/Case